

Precision Dewpoint Meter *with Michell Chilled Mirror Technology*

This reference dewpoint instrument is just one of a wide range of monitoring and calibration instruments offered by Michell. The S8000 combines the ultimate measurement sensitivity and accuracy with flexibility for dynamic measurement and testing protocols or for use as a calibration reference standard.

- Fundamental, accurate and drift free dewpoint measurement
- Measurement range: -75 to +105°F
- $\pm 0.18^\circ\text{F}$ dewpoint accuracy
- Dual optics - supreme sensitivity
- Multiple engineering units
- Data logging to USB or SD card
- Vertical or horizontal configuration



For all your moisture measurement requirements, please see www.michell.com
We will have the solution for you.



Michell Instruments, Inc.
11 Old Sugar Hollow Road
Danbury, CT 06810 USA

Tel: (203) 744-6881 Fax: (203) 743-2051
Email: us.sales@michell.com
Web: www.michell.com